

### 13. Safety classification of ME EQUIPMENT

Protection against electric shock	Internally powered ME equipment
Applied part	Type B applied part, including the whole forehead thermometer
Protection against harmful ingress of water or particulate matter	IPX0
Mode of operation	Continuous operation
<b>Note:</b> Not intended to be sterilized. Not for use in an OXYGEN RICH ENVIRONMENT	

Manufacturer 	Zhongshan JinLi Electronic Weighing Equipment CO.,LTD. 283 <sup>rd</sup> South Min'an Road, Xiaolan Town, Zhongshan City, Guangdong, China Contact Tel.:+86(760)22100980 Fax:+86(760)22108763
EU-representative 	Wellkang Ltd Suite B,29 Harley Street LONDON,W1G 9QR,U.K. www.CE-Marking.EU

# KINLEE

## Non-contact Electronic Forehead Infrared Thermometer FT3010 User Manual



Zhongshan Jinli Electronic Weighing Equipment Co., Ltd.

Ver: 201202A

## Contents

1. Foreword	3
2. Declaration of conformity	3
3. Important Safety Information	3
4. Product Characteristics	6
5. Technical Data	7
6. Intended Use	8
7. Parts Identification (As picture)	9
7.1 Basic components	9
7.2 Display Indicator	10
8. Operation and setting method	10
8.1 Installing or replacing the batteries	10
8.2 Setting method	10
8.3 Advice For The 1 <sup>st</sup> Use	11
8.4 Correct measurement	12
9. Troubleshooting	14
10. Care and Maintenance	14
11. Quality promised and service	15
12. Attachment list	15
13. Safety classification of ME EQUIPMENT	16

- The Infrared Sensor is the most precise part, must be protected carefully.
- Clean and Disinfect the device with a cotton bud lightly moistened with 70% alcohol.
- Do not clean the device with corrosive detergent.
- Keep the device in a dry environment, and keep it away from dust and direct sunlight.
- The main unit is not waterproof. Be careful when handling this unit so that no liquid (alcohol, water, or hot water) will get into the main unit.
- If the infrared sensor becomes dirty, lightly wipe it with a soft dry cloth.
- It is regular recommended to do the recalibration every two years to ensure functioning correct.

**WARNING: Do not modify this equipment without authorization from manufacturer.**

### 11. Quality promised and service

The guarantee time is two years after you buy the device.







**Notice:** If the device is broken by the non-formal used or disassemble by you, the device can not be in the guarantee range.

Tips: Please keep the guarantee card and the buying receipt, which is very useful when you need to amend it.

### 12. Attachment list

- One set of Electronic Forehead Infrared Thermometer.
- One set of user manual.
- One set of the Warranty certificate.
- 2×1.5 AA Batteries(optional).

## 9. Troubleshooting.

Error symbol	Cause	Processing Method
	Low battery voltage	Exchange new batteries
	system fault	Return the selling point
	Ambient temperature too high	Please follow user manual user in the Ambient temperature.
	Ambient temperature too low	
	Measuring temperature too high	Stop measuring
	Measuring temperature too low	

## 10. Care and Maintenance

- The protective glass over the lens is the most important and fragile part of the thermometer, please take great care of it.
- Do not recharge non rechargeable batteries, do not throw in fire.
- Do not expose the thermometer to sunlight or water.

## 1. Foreword


Thank you for purchasing **KINLEC** product. **FT3010 No-contact Electronic Forehead Infrared Thermometer** (here in after called the **Forehead Thermometer** for short) is a thermometer which is used for measuring the temperature of human body at forehead by the principle of receiving infrared. According to the human body skin difference, the measuring temperature would be different. When changing the location from two different temperature location, before use, please put the forehead thermometer in the same condition at least 30 minutes, so that the forehead thermometer could be adapted to the ambient temperature. It is mainly applied to household (consumer use). For your correctly use, please make sure you have read this user manual. Please keep it carefully, in case you can refer again someday.





## 2. Declaration of conformity



- This device fulfils the provisions of EC directive 93/42/EEC (Medical Device Directive).
- This forehead thermometer meets requirements established in ASTM Standard E1965-98(2009). Full responsibility for the conformance of this product to the standard is assumed by Zhongshan Jinli Electronic Weighing Equipment Co., Ltd
- This device has been tested and homologated in accordance with EN60601-1-2:2007 for EMC. This does not guarantee in any way that the device will not be affected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

## 3. Important Safety Information

- The warning symbol and outline mentions in this user manual give guidance of safety and correct method, and avoid harm to your body.
- The warning symbol and outline as follow:

Warning symbol	Meaning
	It may cause harm to body or mangle article.

Cutline	
	Prohibit operation. Prohibit operation can be found in  or description by words or photo nearby. Symbol in left means "General prohibit"
	Must follow operation Must follow operation can be found in  or description by words or photo nearby. Symbol in left means "General must follow operation"

 <b>Batteries use</b>	
<ul style="list-style-type: none"> <li>● Please don't commix new and old batteries or different brand and type of batteries.</li> <li>● Please pull out the batteries if it is not in use for a long time (More than 3 months )</li> <li>● Please set "+、-"pole correctly.</li> <li>● Disobey of above operation may course batteries heat, fluid, or blast, and therefore mangle the product.</li> <li>● If battery fluid should get on your skin or clothing, please rinse with plenty of clean water at once.</li> <li>● Please don't through batteries in fire! It may course blast.</li> <li>● Disposal of discarded batteries should be in accordance with the national regulations for the disposal of discarded batteries.</li> </ul>	

**Cautions:**

- After replaced the batteries or move the thermometer to a new place ,It is suggested to lay the forehead thermometer aside at least 30 minutes before using, so that you can measure temperature more accurately.
- If the measurement is made without an effective distance or an deflection of the central position of the measured object, a deflection from actual temperature may cause. It is suggested to repeat the measurement once or more.
- While measuring, it is suggested to focus the detection window at the forehead, and press the measurement button in front of the handle, scan the forehead around, so that it can measure the highest and lowest temperature of the forehead.
- It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf.

**8.4.3 Selecting measuring body temperature or surface temperature**

Turn on the thermometer, press "SET" key to select body temperature measurement mode or surface temperature measurement mode.

**8.4.4 Data Memory**

In measuring mode, press "+" or "-" key to view measurement result. "+" key for page up, "-" key for page down. Press the "+" or "-" key and hold for 2 seconds, the result switches quickly.

To do it, press the “SET” button for 3 seconds, the screen displays F1, press “SET” button again until you get F3, press “+” button in order to add the difference (in our example, 0.6°C (1.2°F) ).


8.3.3 To check, take the temperature again using the Non-contact body infrared Thermometer.

**Caution: This setting applies only to Body mode.**

## 8.4 Correct measurement


### 8.4.1 Measurement body temperature











Hold the forehead thermometer, leave one finger in front of the button, keep transducer window at human’s forehead, leaving a distance of 50-150mm from the forehead. Press the button, and the thermometer will start up automatically and display the measurement result, as shown in the figure below:









### 8.4.2 Measurement ambient temperature.

Hold the thermometer, leave one finger in front of the button, Level the detection window face to the place where you want to measure, and then press button. The thermometer will start up automatically and display the measurement result, as shown in the figure below:



 <b>Universality affair</b>	
Use of this thermometer is not intended as a substitute for consultation with your physician. There may be a risk of exasperating the state of an illness by self-estimate.	
Please do not put our products close to the charged object. In order to avoid the possibility of electric shock.	
Don't expose this thermometer to extreme temperature conditions of $> 50\text{ }^{\circ}\text{C}$ or $<-20\text{ }^{\circ}\text{C}$ .	
Do not apply a strong shock to, drop, step on, or vibrate the main unit. Do not drop or knock the device, and do not use if damaged.	
Don't use the device beyond intended use.	
Keep the device away from water and heat, including direct sunlight.	
Do not use a portable phone near the unit.	
Do not disassemble, repair, or modify the unit.	
Please use and keep this product with the standard mentions in this manual, Otherwise, It can not measure correctly.	

Meanings of Symbols	
Signs	Notes on the signs
	Type B Applied Part
	Operating instructions
	Manufacturer
	Prohibition of free throw
	European community
	Specific batch code
CE <sub>0120</sub>	CE Mark

#### 4. Product Characteristics

- Designed for measuring the forehead temperature, with a dynamic offset for the ambient temperature and forehead temperature.
- Exclusive use German infrared probe for temperature measurement, with high accuracy and stable performance.
- Providing the function of sound notification of high body temperature.
- With 32 sets of measurement memories.
- LCD digital displayer with backlight.
- Fahrenheit and Celsius optional mode.
- Auto off function to save energy.
- Cute appearance and convenient operation.

#### 8.2.3 Offset value setting---Menu F3 【Default 0.0°C (0.0°F)】

Enter into setting mode, press “SET” key triple, the LCD display “ F3”, then press “+” key to increase the setting value 0.1°C(0.1°F), press “-” key to reduce the setting value 0.1°C(0.1°F).

#### 8.2.4 Sound notification setting

Enter into setting mode, press “SET” key four times, the LCD display “ F4”, then press “+” key to turn on the sound notification function, press “-” key to turn off the sound notification function.

#### 8.2.5 Backlight setting---Menu F5

Enter into setting mode, press “SET” key five times, the LCD display “ F5”, then press “+” key to turn on the backlight function, press “-” key to turn off the backlight function.

#### 8.2.6 Exit setting mode

Press “SET” key until shutdown.

#### 8.3 Advice For The 1<sup>st</sup> Use


For stable and reliable results, it is essential to check the forehead thermometer and change as needed, as follows:

**8.3.1** Take the temperature of the same person using a conventional thermometer; you will get 37.0°C (98.6°F) for instance.

**8.3.2** Take the temperature of the same person using the forehead thermometer keeping the 50 to 150 mm distance between the thermometer and the forehead (Take care to remove any obstacle which could alter the measurement (hair, perspiration...)). If you get 37.0°C (98.6°F), the Non-contact body infrared Thermometer is properly set and ready for use.

If you get a lower temperature, 36.4°C (97.4°F) for example, your difference is 0.6°C (1.2°F), you should adjust the temperature on the forehead Thermometer Body infrared Thermometer and add the difference, i.e. 0.6°C (1.2°F).


## 7.2 Display



## 8. Operation and setting method

### 8.1 Installing or replacing the batteries

a. Push the battery cover back, open the cover from the body.



b. Insert two "AA" batteries as indicated in the battery compartment and then replace the battery cover.

### 8.2 Setting method

Keep pressing "SET" key and hold for 3 seconds into the setting mode.

#### 8.2.1 Select measurement unit---Menu F1

Enter into setting mode, the LCD displays "F1", press "+" key to select "F" unit, press "-" key to select "C" unit.

#### 8.2.2 Warning value setting---Menu F2【Default 37.8°C (100.4°F)】

Enter into setting mode, short press "SET" key twice, the LCD displays " F2", then press "+" key to increase the setting value 0.1°C(0.1°F), press "-" key to reduce the setting value 0.1°C(0.1°F).

## 5. Technical Data

Normal operating conditions	Temperature : 16°C~40°C
	Relative humidity : 15%~95%RH
	860-1060hPa
Transport condition and storage	Power supply : DC 3V (2 x AA LR6)
	Ambient temperature: -20°C ~50°C
	Relative humidity (RH): 15%-95%
Dimension	860-1060hPa
	90mm*43mm*148mm(length*width*height)
Weight (net weight)	120g (not including the batteries)
LCD display resolution	0.1°C (0.1°F)
Displayed Temperature Range	Body temperature:22.0~42.9°C (71.6~109.0°F)
	Surface temperature : 0~100°C(32.0~212.0°F)
Maximum Laboratory Error	1) 22.0~42.0 °C : 0.3°C
	2) out of 22.0~42.0 °C :0.5°C
Power	≤50mW
Measurement time	≤ 1second
Measurement distance	50mm~150mm
Blackbody recommended for verifying	The blackbody addressed in 6.1.3.3 of E1965-98(2009)
Auto-off time	10 second
measurement data memories	32 sets
Body temperature Mode(adjusted)	The displayed value display body temperature adjusted to auxiliary.
Surface temperature Mode(not adjusted)	The displayed value display surface temperature (Forehead temperature or object surface temperature)

**WARNING:** the performance of the instrument may be adversely affected the following occur:

- Operation outside of the specified subject temperature t.
- Operation outside of the specified ambient temperature and humidity ranges
- Storage outside of the specified ambient temperature and humidity ranges
- Mechanical shock may be adversely affected performance of the instrument.

## 6. Intended Use

The forehead thermometer is designed for body surface and forehead temperature measurement for infant and adults without contact to human body.

The forehead thermometer can also be used to measure the temperature of a baby-bottle or bath or room temperature (by using the surface mode).

Normal Temperatures According To Measurement Method

Measurement Method	Normal Temp °C	Normal Temp °F
Rectal	36.6 to 38	97.8 to 100.4
Oral	35.5 to 37.5	95.9 to 99.5
Axillary	34.7 to 37.3	94.4 to 99.1
Ear	35.8 to 38	96.4 to 100.4

The temperature of the human body varies throughout the day. It can also be influenced by numerous external factors: age, sex, type and thickness of skin...

Normal Temperatures According To Age

Age	Temp °C	Temp °F
0-2 years	36.4 to 38.0	97.5 to 100.4
3-10 years	36.1 to 37.8	97.0 to 100.0
11-65 years	35.9 to 37.6	96.6 to 99.7
> 65 years	35.8 to 37.5	96.4 to 99.5

**Attention:** ASTM laboratory accuracy requirements in the display of 37 to 39°C (98 to 102°F) for Skin IR thermometers is  $\pm 0.3^{\circ}\text{C}$  ( $\pm 0.5^{\circ}\text{F}$ ), whereas for mercury in-glass and electronic thermometers, the requirement per ASTM Standards E667-86 and 1112-86 is  $\pm 0.1^{\circ}\text{C}$  ( $\pm 0.2^{\circ}\text{F}$ ).

**When you need to clinical accuracy characteristics and procedures, please contact with the local distributors obtain from our company.**

## 7. Parts Identification (As picture)

### 7.1 Basic components

